Codes: 202811, 202810, 202808, 202809, 202822, 202801, 202807, 202802, 202814, 202803, 202821, 202813, 202818, 202804, 202815, 202820, 202839
Polyethylene wide mouth jars with lid and insert plug

2006 catalogue: page 142, 250, 251

Code	Capacity ml	Unitary weight without cap g	Internal mouth Φ mm	Body Φ mm	Height with cap mm	Case Quantity	Case weight	Case volume
202811	30	4,2	25	32	53	100	1.00	0.011
202810	60	8,2	38	52	48	200	4.00	0.049
202808	90	8	51	60	51	150	3.80	0.049
202809	125	12	51	60	62	200	5.00	0.087
202822	135	11,5	40	53	81	180	5.00	0.087
202801	170	14	51	60	80	160	4.70	0.087
202807	200	19,7	49	64	94	150	5.00	0.087
202802	250	20	54	67	100	125	4.70	0.087
202814	400	27,3	60	74	124	130	9.00	0.140
202803	500	35	67	80	131	120	8.40	0.140
202821	500	44,9	86	103	93	95	9.00	0.140
202813	750	48	86	103	127	75	7.70	0.140
202818	1.000	56,7	86	103	157	50	5.80	0.140
202804	1.250	56,7	86	103	173	50	6.10	0.140
202815	1.550	64,6	86	113	181	45	6.00	0.140
202820	2.000	83,1	86	103	283	35	5.80	0.140
202839	2.000	73	105	128	205	32	4.00	0.140

Medical Device IVDD 98/79/CE

SPECIFIC USE AND DESCRIPTION

Polyethylene jars with screw cap and insert plug Jars and plugs are natural translucent colour; lids are black. Ask for caps and jars in other colour.


Supplied uncapped (excepting code 202839 which is supplied capped)

The material (high density polyethylene) used for the production of this article fulfills the requirements to come into contact with food.

Upon request, jars can be supplied individually wrapped for a minimum order of 3 cases, and sterile for a minimum order of 6 cases.

Elaborated by Marketing	Verified by Quality	Authorized by Management





Elaborated by Marketing	Verified by Quality	Authorized by Management