
3050 Spruce Street, Saint Louis, MO 63103, USA
Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
Eosin Y disodium salt - certified by the Biological Stain Commission

Product Number:	**E4382**
CAS Number:	17372-87-1
MDL:	MFCD00005040
Formula:	C20H6Br4Na2O5
Formula Weight:	691.85 g/mol



TEST	Specification
Appearance (Form)	Powder
Infrared Spectrum	Conforms to Structure
Solubility (Color)	Very Light Orange to Dark Orange to Dark Red
c = 1mg/ml; Water	
Solubility (Turbidity)	Clear
Biological Stain Commission	Pass
Commission Certified Number:	
Dye Content	≥ 90 %

Specification: PRD.4.ZQ5.10000042329

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.