
Technical Drawing Pharma Focus Tubing (TPF)

FREUDENBERG MEDICAL

Liveo C6- 150

Materialnumber	OD (B)	±	ID (A)	±	Wall (C)	Length
00 45 63 25 21- (61-252-21)	6,35 mm	0,18 mm	3,18 mm	0,18 mm	1,59 mm	15 m
00 45 63 25 22- (61-252-22)	7,95 mm	0,25 mm	4,78 mm	0,25 mm	1,59 mm	15 m
00 45 63 25 23- (61-252-23)	9,53 mm	0,25 mm	4,78 mm	0,25 mm	2,38 mm	15 m
00 45 63 25 25- (61-252-25)	11,13 mm	0,38 mm	6,35 mm	0,25 mm	2,39 mm	15 m
00 45 63 25 26- (61-252-26)	12,71 mm	0,38 mm	6,35 mm	0,25 mm	3,18 mm	15 m
00 45 63 25 27- (61-252-27)	12,71 mm	0,38 mm	7,95 mm	0,25 mm	2,38 mm	15 m
00 45 63 25 30- (61-252-30)	15,88 mm	0,51 mm	9,53 mm	0,38 mm	3,18 mm	15 m
00 45 63 25 32- (61-252-32)	19,06 mm	0,51 mm	12,71 mm	0,38 mm	3,18 mm	15 m
00 45 63 25 33- (61-252-33)	22,23 mm	0,51 mm	15,88 mm	0,51 mm	3,18 mm	15 m
00 45 63 25 35- (61-252-35)	25,41 mm	0,51 mm	19,06 mm	0,51 mm	3,18 mm	15 m
00 45 63 25 36- (61-252-36)	38,12 mm	0,76 mm	25,41 mm	0,51 mm	6,36 mm	15 m



Kaiserslautern, 08.12.2021

Quality Control